



---

## Product Datasheet

### CXXC4 Antibody (orb1249423)

<b>Description</b>	CXXC4 Antibody	 <p>250kDa 150kDa 100kDa</p>
<b>Species/Host</b>	Goat	 <p>50kDa 37kDa 25kDa 20kDa 15kDa</p> <p>orb1249423 (3.8 ug/ml) staining of paraf...</p>
<b>Reactivity</b>	Human	
<b>Conjugation</b>	Unconjugated	
<b>Tested Applications</b>	ELISA, IHC	
<b>Immunogen</b>	The immunogen for this antibody is: C-GTSLERTPVPSAEA	
<b>Target</b>	CXXC4	
<b>Preservatives</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.	
<b>Form/Appearance</b>	Liquid	
<b>Concentration</b>	500 ug/mL	
<b>Storage</b>	Aliquot and store at -20°C. Minimize freezing and thawing.	
<b>Note</b>	For research use only	
<b>Application notes</b>	<p>Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Preliminary experiments in Human Brain (Cerebellum, Hippocampus) lysates gave no specific signal but low background (at antibody concentration up to 1ug/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Immunohistochemistry: In paraffin embedded Human Breast shows cytoplasm staining in myoepithelial cells. Recommended concentration, 3-5ug/ml.</p>	
<b>Clonality</b>	Polyclonal	
<b>Uniprot ID</b>	<a href="#">Q9H2H0</a>	
<b>NCBI</b>	<a href="#">NP_079488.1</a>	
<b>Dilution Range</b>	<p>Peptide ELISA: antibody detection limit dilution 1:32000. Western Blot: Preliminary experiments in Human Brain (Cerebellum, Hippocampus) lysates gave no specific signal but low background (at antibody concentration up to 1ug/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Immunohistochemistry: In paraffin embedded Human Breast shows cytoplasm staining in myoepithelial cells. Recommended concentration, 3-5ug/ml.</p>	
<b>Expiration Date</b>	12 months from date of receipt.	