
Product Datasheet

TMPRSS4 Antibody (orb1249595)

Description

Tmprss4 Antibody

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications
Immunogen

The immunogen for this antibody is: HCFRKHTDVFNWK

Target

Tmprss4

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Form/Appearance

Liquid

Concentration

500 ug/mL

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

Peptide ELISA: antibody detection limit dilution 1:1000. Western Blot: Approx. 48kDa band observed in lysates of cell lines HEK293 and A549 (calculated MW of 48.2kDa according to NP_063947.1). Recommended concentration: 0.3-1ug/ml.

Clonality


Polyclonal

Uniprot ID
[Q9NRS4](#)
NCBI
[NP_001077416.1](#), [NP_063947.1](#)
Dilution Range

Peptide ELISA: antibody detection limit dilution 1:1000. Western Blot: Approx. 48kDa band observed in lysates of cell lines HEK293 and A549 (calculated MW of 48.2kDa according to NP_063947.1). Recommended concentration: 0.3-1ug/ml.

Expiration Date

12 months from date of receipt.


 orb1249595
(1 ug/ml)
staining of
HEK293 ...