


Product Datasheet

HAMP Antibody (orb1249635)

Description	HAMP Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	HAMP (AAH20612, 25 a.a. ~ 85 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Target	HAMP
Preservatives	PBS, pH 7.2
Form/Appearance	Liquid
Concentration	0.48 mg/ml
Storage	Store Hepcidin antibody at 4 °C or -20 °C. As with all antibodies avoid freeze/thaw cycles.
Note	For research use only
Application notes	Hepcidin antibody can be used in Western Blot, immunohistochemistry starting at 1:50, and immunoprecipitation.
Isotype	IgG1,k
Clonality	Monoclonal
Uniprot ID	P81172
NCBI	P81172
Dilution Range	Hepcidin antibody can be used in Western Blot, immunohistochemistry starting at 1:50, and immunoprecipitation.
Expiration Date	12 months from date of receipt.



Immunohistochemistry of human liver tiss...