


## Product Datasheet

### STK17A Antibody (orb1249714)

<b>Description</b>	STK17A Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC, RNA-I, WB
<b>Immunogen</b>	STK17A (AAH47696, 301 a.a. ~ 414 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
<b>Target</b>	STK17A
<b>Preservatives</b>	PBS, pH 7.2
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	0.3 mg/ml
<b>Storage</b>	Store DRAK1 antibody at -20 °C. Aliquot to avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	DRAK1 antibody can be used in Western Blot, and immunohistochemistry starting at 5 µg/mL.
<b>Isotype</b>	IgG2a,k
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">Q9UEE5</a>
<b>NCBI</b>	<a href="#">Q9UEE5</a>
<b>Dilution Range</b>	DRAK1 antibody can be used in Western Blot, and immunohistochemistry starting at 5 µg/mL.
<b>Expiration Date</b>	12 months from date of receipt.



Immunohistochemistry of human placenta t...