


Product Datasheet

SPANXC Antibody (orb1249762)

Description	SPANXC Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	SPANXC antibody was raised against amino acids from the C-Terminus of SPANXC (Human) Amino acids 31-47.
Target	SPANXC
Preservatives	0.02 M potassium phosphate, 0.15 M sodium chloride, pH 7.2, 0.01% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Short term: 4°C. Long term: Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	SPANXC antibody can be used in Western Blot, and immunohistochemistry starting at 10 µg/mL.
Clonality	Polyclonal
Uniprot ID	Q9NY87
NCBI	Q9NY87
Dilution Range	SPANXC antibody can be used in Western Blot, and immunohistochemistry starting at 10 µg/mL.
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of SPANXC ...