
Product Datasheet

PPP1R13L Antibody (orb1249766)

Description

PPP1R13L Antibody


Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IF, IHC, WB

Applications

Immunohistochemistry staining of PPP1R13...

Immunogen

Amino acids 780-797 of human iASPP (isoform 1) protein.

Target

PPP1R13L

Preservatives

0.02 M potassium phosphate, 0.15 M sodium chloride, pH 7.2, 0.01% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/ml

Storage

Store PPP1R13L antibody at -20 °C prior to opening. Dilute only prior to immediate use. For extended storage aliquot contents and freeze at -20 °C or below. As with all antibodies avoid freeze/thaw cycles.

Note

For research use only

Application notes

PPP1R13L antibody can be used in Western Blot starting at 1 - 3 µg/mL, immunohistochemistry starting at 5 µg/mL, and flow cytometry.

Clonality

Polyclonal

Uniprot ID
[Q8WUF5](#)
NCBI
[Q8WUF5](#)
Dilution Range

PPP1R13L antibody can be used in Western Blot starting at 1 - 3 µg/mL, immunohistochemistry starting at 5 µg/mL, and flow cytometry.

Expiration Date

12 months from date of receipt.