
Product Datasheet

EGR1 Antibody (orb1249781)

Description

EGR1 Antibody


Species/Host

Rabbit

Reactivity

Human, Mouse, Primate

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunohistochemistry staining of Early G...

Immunogen

Early Growth Response 1 antibody was raised against amino acids 94 - 108 of Early Growth Response 1 (Human).

Target

EGR1

Preservatives

0.02 M potassium phosphate, 0.15 M sodium chloride, pH 7.2, 0.01% sodium azide.

Form/Appearance

Liquid

Concentration

0.93 mg/ml

Storage

Store Early Growth Response 1 antibody at 4 °C or -20 °C. As with all antibodies avoid freeze/thaw cycles.

Note

For research use only

Application notes

Early Growth Response 1 antibody can be used in ELISA, Western Blot starting at 2 - 4 µg/mL, immunohistochemistry starting at 5 µg/mL, and immunofluorescence.

Clonality

Polyclonal

Uniprot ID
[P18146](#)
NCBI
[P18146](#)
Dilution Range

Early Growth Response 1 antibody can be used in ELISA, Western Blot starting at 2 - 4 µg/mL, immunohistochemistry starting at 5 µg/mL, and immunofluorescence.

Expiration Date

12 months from date of receipt.