
Product Datasheet

MLF1IP Antibody (orb1249782)

Description

MLF1IP Antibody

Species/Host

Rabbit

Reactivity


Bovine, Canine, Human, Primate

Conjugation

Unconjugated

Tested

ELISA, IHC, IP, WB

Applications


Immunohistochemistry staining of MLF1IP ...

Immunogen

MLF1IP antibody was raised against amino acids surrounding Thr78 of MLF1IP (Human).

Target

MLF1IP

Preservatives

0.02 M potassium phosphate, 0.15 M sodium chloride, pH 7.2, 0.01% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/ml

Storage

Store MLF1IP antibody at -20 °C prior to opening. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20 °C or below. As with all antibodies avoid freeze/thaw cycles.

Note

For research use only

Application notes

MLF1immunoprecipitation antibody can be used in ELISA, Western Blot, immunohistochemistry starting at 5 µg/mL, and immunoprecipitation

Clonality

Polyclonal

Uniprot ID
[Q71F23](#)
NCBI
[Q71F23](#)
Dilution Range

MLF1immunoprecipitation antibody can be used in ELISA, Western Blot, immunohistochemistry starting at 5 µg/mL, and immunoprecipitation

Expiration Date

12 months from date of receipt.