


Product Datasheet

NCS1 Antibody (orb1249784)

Description	NCS1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Frequenin antibody was raised against recombinant protein from mouse with extensive post-translational modifications (Mouse).
Target	NCS1
Preservatives	0.02 M potassium phosphate, 0.15 M sodium chloride, pH 7.2, 0.01% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store Frequenin antibody at -20 °C prior to opening. Dilute only prior to immediate use. For extended storage aliquot contents and freeze at -20 °C or below. As with all antibodies avoid freeze/thaw cycles.
Note	For research use only
Application notes	Frequenin antibody can be used in ELISA, and immunohistochemistry starting at 20 µg/mL.
Clonality	Polyclonal
Uniprot ID	P62166
NCBI	P62166
Dilution Range	Frequenin antibody can be used in ELISA, and immunohistochemistry starting at 20 µg/mL.
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of Frequen...