

## Product Datasheet

### CCL20 Antibody (orb1249863)

<b>Description</b>	CCL20 Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Immunogen</b>	E.coli derived recombinant Human MIP-3alpha.
<b>Target</b>	CCL20
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store Lyophilized at room temperature up to 1 month; Reconstituted for up to two weeks at 2-8°C. Aliquot and freeze at -20°C for long term storage. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	CCL20 antibody can be used for ELISA starting at 2 - 4 µg/mL, staining paraffin embedded tissues starting at 10 µg/mL, and Western Blot starting at 0.2 - 0.4 µg/mL.
<b>Isotype</b>	IgG1,k
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P78556</a>
<b>NCBI</b>	<a href="#">P78556</a>
<b>Dilution Range</b>	CCL20 antibody can be used for ELISA starting at 2 - 4 µg/mL, staining paraffin embedded tissues starting at 10 µg/mL, and Western Blot starting at 0.2 - 0.4 µg/mL.
<b>Expiration Date</b>	12 months from date of receipt.