


Product Datasheet

MKRN1 antibody (orb125005)

Description	Goat polyclonal antibody to MKRN1
Species/Host	Goat
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Makorin
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: In transfected HEK293 transiently expressing full-length Human MKRN1 (myc and DYKDDDDK tagged), a band of approx. 60kDa was observed. No bands were observed in mock-transfected HEK293 and the same band was observed using anti-DYKDDDDK tag antibody. Recom
Protein Sequence	VTKEESEKEQTAVET
Clonality	Polyclonal
MW	53.4; 35.2
Entrez	23608
Dilution Range	ELISA: 1:8000, WB: 1-2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 lysate u...