


## Product Datasheet

### HCRT Antibody (orb1250096)

<b>Description</b>	HCRT Antibody
<b>Species/Host</b>	Sheep
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC
<b>Immunogen</b>	Prepro - orexin antibody was raised against a synthetic peptide corresponding to C-Terminus of Prepro - orexin (Human).
<b>Target</b>	HCRT
<b>Preservatives</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	9.7 mg/ml
<b>Storage</b>	Store Prepro-orexin antibody at 4 °C or -20 °C. As with all antibodies avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Prepro-orexin antibody can be used in Western Blot, immunohistochemistry starting at 10 - 20 µg/mL, and immunoprecipitation.
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">O43612</a>
<b>NCBI</b>	<a href="#">O43612</a>
<b>Dilution Range</b>	Prepro-orexin antibody can be used in Western Blot, immunohistochemistry starting at 10 - 20 µg/mL, and immunoprecipitation.
<b>Expiration Date</b>	12 months from date of receipt.



Immunohistochemistry staining of Prepro-...