


Product Datasheet

GOT1 Antibody (orb1250097)

Description	GOT1 Antibody
Species/Host	Sheep
Reactivity	Human, Porcine
Conjugation	Unconjugated
Tested Applications	IF, IHC, IP, WB
Immunogen	GOT2 antibody was raised against Aspartate Transaminase from pig heart (Porcine).
Target	GOT1
Preservatives	0.015 M potassium phosphate buffer, 0.15 M sodium chloride, pH 7.2, 0.1% sodium azide.
Form/Appearance	Liquid
Concentration	10 mg/ml
Storage	Short term: -20°C; Long term: -20°C.
Note	For research use only
Application notes	GOT2 antibody can be used in Western Blot, immunohistochemistry starting at 5 µg/mL, and immunoprecipitation.
Clonality	Polyclonal
Uniprot ID	P17174
NCBI	P17174
Dilution Range	GOT2 antibody can be used in Western Blot, immunohistochemistry starting at 5 µg/mL, and immunoprecipitation.
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of GOT2 in...