

## Product Datasheet

### SLC39A14 Antibody (orb1250145)

<b>Description</b>	SLC39A14 Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Gorilla, Human, Monkey, Porcine, Primate
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC
<b>Immunogen</b>	Synthetic 15 amino acid peptide from cytoplasmic domain of human SLC39A14.
<b>Target</b>	SLC39A14
<b>Preservatives</b>	PBS, 0.1% sodium azide.
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	SLC39A14 Antibody can be used for Immunohistochemistry in paraffin embedded tissues starting at 5 - 7.5 µg/mL.
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q15043</a>
<b>NCBI</b>	<a href="#">Q15043</a>
<b>Dilution Range</b>	SLC39A14 Antibody can be used for Immunohistochemistry in paraffin embedded tissues starting at 5 - 7.5 µg/mL.
<b>Expiration Date</b>	12 months from date of receipt.