


Product Datasheet

PAG1 Antibody (orb1250273)

Description	PAG1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	PAG1 antibody was raised against amino acids from the C-Terminus of PAG1 (Human).
Target	PAG1
Preservatives	PBS, 15 mM sodium azide, approx., pH 7.4.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store PAG1 antibody at 2-8 °C. Do not freeze.
Note	For research use only
Application notes	IHC-P (10 ug/ml), IP, WB
Clonality	Polyclonal
Uniprot ID	Q9NWQ8
NCBI	Q9NWQ8
Dilution Range	IHC-P (10 ug/ml), IP, WB
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of PAG1 in...