


Product Datasheet

CSF1R antibody (orb125029)

Description	Goat polyclonal antibody to CSF1R
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CSF1R (aa734-744)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: In transfected HEK293 transiently expressing full-length Human CSF1R (myc and DYKDDDDK tagged), a band of approx. 200kDa was observed. No bands were observed in mock-transfected HEK293 and the same band was observed using anti-myc tag antibody. Recommend
Protein Sequence	NDSFSEQDLDK
Clonality	Polyclonal
MW	108
Entrez	1436
Dilution Range	ELISA: 1:8000, WB: 0.2-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 lysate u...