


Product Datasheet

DCN Antibody (orb1250309)

Description	DCN Antibody
Species/Host	Sheep
Reactivity	Bovine, Gallus, Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Decorin antibody was raised against a synthetic peptide corresponding to C-Terminus of Decorin (Human).
Target	DCN
Preservatives	PBS, pH 7.2, 0.1% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Store Decorin antibody at 4 °C or -20 °C. As with all antibodies avoid freeze/thaw cycles.
Note	For research use only
Application notes	Decorin antibody can be used in immunohistochemistry starting at 10 µg/mL.
Clonality	Polyclonal
Uniprot ID	P07585
NCBI	P07585
Dilution Range	Decorin antibody can be used in immunohistochemistry starting at 10 µg/mL.
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of DCN in ...