

Product Datasheet

GRAMD3 antibody (orb125040)

Description	Goat Polyclonal to GRAMD3. NP_001139791.1; NP_076416.2; NP_001139792.1; NP_001139793.1; NP_001139794...
Species/Host	Goat
Reactivity	Bovine, Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	GRAMD3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments in Human Brain (frontal Cortex), Spleen and Placenta lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml).
Protein Sequence	EAKPTRADAH
Clonality	Polyclonal
MW	49.9; 47.9; 36.6; 49.6; 46.3
Entrez	65983
Dilution Range	ELISA: 1:1000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.