


Product Datasheet

PDE1C Antibody (orb1250452)

Description	PDE1C Antibody
Species/Host	Rabbit
Reactivity	Bat, Bovine, Canine, Equine, Gorilla, Hamster, Human, Monkey, Mouse, Primate, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	PDE1C antibody was raised against a peptide located in the C-Terminal domain of PDE1C (Human).
Target	PDE1C
Preservatives	PBS, 0.1% sodium azide.
Form/Appearance	Liquid
Concentration	0.21 mg/ml
Storage	PDE1C antibody should be stored long term (months) at -80 °C and short term (days) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
Note	For research use only
Application notes	PDE1C antibody can be used in immunohistochemistry starting at 10 µg/mL, and flow cytometry.
Clonality	Polyclonal
Uniprot ID	Q14123
NCBI	Q14123
Dilution Range	PDE1C antibody can be used in immunohistochemistry starting at 10 µg/mL, and flow cytometry.
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of PDE1C i...