


Product Datasheet

FLT3 Antibody (orb1250458)

Description	FLT3 Antibody
Species/Host	Rabbit
Reactivity	Gorilla, Human, Primate
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	FLT3 antibody was raised against a peptide located near the internal domain of FLT3 (Human).
Target	FLT3
Preservatives	PBS, 0.1% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	FLT3 antibody should be stored long term (months) at -80 °C and short term (days) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
Note	For research use only
Application notes	FLT3 antibody can be used in immunohistochemistry starting at 13 - 26 µg/mL.
Clonality	Polyclonal
Uniprot ID	P36888
NCBI	P36888
Dilution Range	FLT3 antibody can be used in immunohistochemistry starting at 13 - 26 µg/mL.
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of FLT3 in...