


Product Datasheet

GPR84 Antibody (orb1250498)

Description	GPR84 Antibody
Species/Host	Rabbit
Reactivity	Gorilla, Human, Primate
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	GPR84 antibody was raised against a peptide located in the 3rd cytoplasmic domain of GPR84 (Human).
Target	GPR84
Preservatives	PBS, 0.1% sodium azide.
Form/Appearance	Liquid
Concentration	0.78 mg/ml
Storage	GPR84 antibody should be stored long term (months) at -80 °C and short term (days) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
Note	For research use only
Application notes	ELISA IHC-P (5 - 10 ug/ml)
Clonality	Polyclonal
Uniprot ID	Q9NQS5
NCBI	Q9NQS5
Dilution Range	ELISA IHC-P (5 - 10 ug/ml)
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of GPR84 i...