


Product Datasheet

FCGR3A Antibody (Biotin) (orb1250664)

Description	FCGR3A Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC
Immunogen	Mononuclear cells from a prolymphocytic B-leukemia
Target	FCGR3A
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD16 Antibody [GRM1] for use in flow cytometry, immunohistochemistry, western blotting, immunoprecipitation, and ELISA assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P08637
NCBI	P08637
Dilution Range	CD16 Antibody [GRM1] for use in flow cytometry, immunohistochemistry, western blotting, immunoprecipitation, and ELISA assays.
Expiration Date	12 months from date of receipt.



Human granulocytes were stained with Mou...