


Product Datasheet

ITGAX Antibody (PE) (orb1250684)

Description	ITGAX Antibody (PE)
Species/Host	Mouse
Reactivity	Human, Primate
Conjugation	PE
Tested Applications	FC
Immunogen	Rheumatoid synovial cells and human monocytes
Target	ITGAX
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD11c Antibody [3.9] for use in flow cytometry assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P20702
NCBI	P20702
Dilution Range	CD11c Antibody [3.9] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



Human lymphocytes, monocytes, and granul...