


Product Datasheet

ITGAM Antibody (Biotin) (orb1250690)

Description	ITGAM Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human, Porcine, Primate
Conjugation	Biotin
Tested Applications	Blocking, Cell Separation, FC, ICC, IP, WB
Immunogen	Rheumatoid synovial cells and human monocytes
Target	ITGAM
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD11b Antibody [ICRF44] for use in flow cytometry, immunohistochemistry / immunocytochemistry, western blotting, and immunoprecipitation assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P11215
NCBI	P11215
Dilution Range	CD11b Antibody [ICRF44] for use in flow cytometry, immunohistochemistry / immunocytochemistry, western blotting, and immunoprecipitation assays.
Expiration Date	12 months from date of receipt.



Human peripheral blood granulocytes were...