

## Product Datasheet

### MS4A1 Antibody (Biotin) (orb1250768)

<b>Description</b>	MS4A1 Antibody (Biotin)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Feline, Human
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	FC
<b>Target</b>	MS4A1
<b>Preservatives</b>	Supplied in PBS/NaN <sub>3</sub>
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store vial at 2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	CD20 Antibody [B-Ly1] for use in flow cytometry, immunohistochemistry, immunoprecipitation, and separation assays.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P11836</a>
<b>NCBI</b>	<a href="#">P11836</a>
<b>Dilution Range</b>	CD20 Antibody [B-Ly1] for use in flow cytometry, immunohistochemistry, immunoprecipitation, and separation assays.
<b>Expiration Date</b>	12 months from date of receipt.