

Product Datasheet

MS4A1 Antibody (FITC) (orb1250769)

Description	MS4A1 Antibody (FITC)
Species/Host	Mouse
Reactivity	Feline, Human
Conjugation	FITC
Tested Applications	FC
Target	MS4A1
Preservatives	Supplied in PBS/NaN ₃
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD20 Antibody [B-Ly1] for use in flow cytometry and immunohistochemistry assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P11836
NCBI	P11836
Dilution Range	CD20 Antibody [B-Ly1] for use in flow cytometry and immunohistochemistry assays.
Expiration Date	12 months from date of receipt.