


Product Datasheet

CD19 Antibody (APC/Cy7) (orb1250772)

Description	CD19 Antibody (APC/Cy7)
Species/Host	Mouse
Reactivity	Human
Conjugation	APC/Cy7
Tested Applications	FC, IHC, IP
Immunogen	NALM-1 and NALM-16 leukemia cell line
Target	CD19
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD19 Antibody [SJ25-C1] for use in flow cytometry assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P15391
NCBI	P15391
Dilution Range	CD19 Antibody [SJ25-C1] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



Human peripheral blood lymphocytes were ...