

Product Datasheet

CD19 Antibody (APC/Cy5.5) (orb1250773)

Description	CD19 Antibody (APC/Cy5.5)
Species/Host	Mouse
Reactivity	Human
Conjugation	APC/Cy5.5
Tested Applications	FC
Immunogen	NALM-1 and NALM-16 leukemia cell line
Target	CD19
Preservatives	Supplied in PBS/NaN ₃ and a stabilizing agent
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD19 Antibody [SJ25-C1] for use in flow cytometry assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P15391
NCBI	P15391
Dilution Range	CD19 Antibody [SJ25-C1] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.