


Product Datasheet

NP Antibody (HRP) (orb1250803)

Description	NP Antibody (HRP)
Species/Host	Mouse
Reactivity	Virus
Conjugation	HRP
Tested Applications	Activity Assay, FC, IHC-Fr, IP
Immunogen	Recombinant influenza virus type B nucleoprotein
Target	NP
Preservatives	Supplied in a stock solution of 50% glycerol/50% PBS, pH 7.4; No preservative added
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	Influenza B, Nucleoprotein Antibody [1H7a] for use in ELISA and western blotting assays.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Influenza B, Nucleoprotein Antibody [1H7a] for use in ELISA and western blotting assays.
Expiration Date	12 months from date of receipt.



Recombinant influenza proteins were reso...