

## Product Datasheet

# Influenza A, Nucleoprotein Antibody [FluA-NP 4F1] (HRP) (orb1250805)

<b>Description</b>	Influenza A, Nucleoprotein Antibody [FluA-NP 4F1] (HRP)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	FC
<b>Immunogen</b>	Recombinant influenza virus type A nucleoprotein
<b>Preservatives</b>	Supplied in a stock solution of 50% glycerol/50% PBS, pH 7.4; No preservative added
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store vial at 2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Influenza A, Nucleoprotein Antibody [FluA-NP 4F1] for use in immunocytochemistry, ELISA, and western blotting assays.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	Influenza A, Nucleoprotein Antibody [FluA-NP 4F1] for use in immunocytochemistry, ELISA, and western blotting assays.
<b>Expiration Date</b>	12 months from date of receipt.