

## Product Datasheet

### S-Tag Antibody [SBSTAGa] (HRP) (orb1250813)

<b>Description</b>	S-Tag Antibody [SBSTAGa] (HRP)
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	FC, IHC-Fr
<b>Immunogen</b>	Recombinant S-tag protein
<b>Preservatives</b>	Supplied in a stock solution of 50% glycerol/50% PBS, pH 7.4; No preservative added
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store vial at 2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	S-Tag Antibody [SBSTAGa] for use in ELISA and western blotting assays.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	S-Tag Antibody [SBSTAGa] for use in ELISA and western blotting assays.
<b>Expiration Date</b>	12 months from date of receipt.