


Product Datasheet

IL33 Antibody (Biotin) (orb1250826)

Description	IL33 Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	AM, ELISA, FC, ICC, IHC, IHC-P, IP, WB
Immunogen	E. coli-expressed human IL-33
Target	IL33
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	IL-33 Antibody [SB127c] for use in ELISA assays.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	095760
NCBI	095760
Dilution Range	IL-33 Antibody [SB127c] for use in ELISA assays.
Expiration Date	12 months from date of receipt.



Standard curve generated with Mouse Anti...