


## Product Datasheet

### MAPK1, MAPK3 Antibody (orb1250831)

<b>Description</b>	MAPK1, MAPK3 Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FLISA, IHC-Fr, IHC-P
<b>Immunogen</b>	Recombinant human Erk1 and Erk2
<b>Target</b>	MAPK1, MAPK3
<b>Preservatives</b>	Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	0.5 mg/mL
<b>Storage</b>	Store vial at 2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	ERK1/2 Antibody [SB46b] for use in ELISA and western blotting assays.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P27361</a> , <a href="#">P28482</a>
<b>NCBI</b>	<a href="#">P28482</a> , <a href="#">P27361</a>
<b>Dilution Range</b>	ERK1/2 Antibody [SB46b] for use in ELISA and western blotting assays.
<b>Expiration Date</b>	12 months from date of receipt.



Total cell lysates from Jurkat cells wer...