

## Product Datasheet

### NP Antibody (orb1250839)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | NP Antibody   |
| <b>Species/Host</b>        | Mouse   |
| <b>Reactivity</b>          | Virus   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | FC, IHC-Fr  |
| <b>Immunogen</b>           | Recombinant influenza virus type B nucleoprotein  |
| <b>Target</b>              | NP  |
| <b>Preservatives</b>       | Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Concentration</b>       | 0.5 mg/mL   |
| <b>Storage</b>             | Store vial at 2-8°C   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Influenza B, Nucleoprotein Antibody [1H7a] for use in ELISA and western blotting assays.            |
| <b>Isotype</b>             | IgG1  |
| <b>Clonality</b>           | Monoclonal  |
| <b>Dilution Range</b>      | Influenza B, Nucleoprotein Antibody [1H7a] for use in ELISA and western blotting assays.            |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |