

Product Datasheet

NP Antibody (orb1250840)

Description	NP Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Recombinant influenza virus type B nucleoprotein
Target	NP
Preservatives	Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	Influenza B, Nucleoprotein Antibody [3D7] for use in ELISA and western blotting assays.
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	Influenza B, Nucleoprotein Antibody [3D7] for use in ELISA and western blotting assays.
Expiration Date	12 months from date of receipt.