


Product Datasheet

Influenza A, Nucleoprotein Antibody [FluA-NP 2C9] (orb1250842)

Description	Influenza A, Nucleoprotein Antibody [FluA-NP 2C9]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	Activity Assay, FC, IHC-Fr, IP
Immunogen	Recombinant influenza virus type A nucleoprotein
Preservatives	Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	Influenza A, Nucleoprotein Antibody [FluA-NP 2C9] for use in immunohistochemistry, ELISA, and western blotting assays.
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	Influenza A, Nucleoprotein Antibody [FluA-NP 2C9] for use in immunohistochemistry, ELISA, and western blotting assays.
Expiration Date	12 months from date of receipt.



Recombinant influenza proteins were reso...