


Product Datasheet

IL3 Antibody (Biotin) (orb1250860)

Description	IL3 Antibody (Biotin)
Species/Host	Rat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC
Immunogen	COS-expressed mouse IL-3
Target	IL3
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	IL-3 Antibody [MP2-43D11] for use in immunohistochemistry / immunocytochemistry, western blotting, ELISA, ELISpot, and multiplex assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P01586
NCBI	P01586
Dilution Range	IL-3 Antibody [MP2-43D11] for use in immunohistochemistry / immunocytochemistry, western blotting, ELISA, ELISpot, and multiplex assays.
Expiration Date	12 months from date of receipt.



Standard curve generated with Rat Anti-M...