


## Product Datasheet

### IL4 Antibody (Biotin) (orb1250878)

<b>Description</b>	IL4 Antibody (Biotin)
<b>Species/Host</b>	Rat
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	FC
<b>Immunogen</b>	E. coli-expressed mouse IL-4 (glutaraldehyde crosslinked)
<b>Target</b>	IL4
<b>Preservatives</b>	Supplied in PBS/NaN3
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	0.5 mg/mL
<b>Storage</b>	Store vial at 2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	IL-4 Antibody [BVD6-24G2] for use in flow cytometry, immunohistochemistry, western blotting, immunoprecipitation, ELISA, ELISpot, and multiplex assays.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P07750</b>
<b>NCBI</b>	<b>P07750</b>
<b>Dilution Range</b>	IL-4 Antibody [BVD6-24G2] for use in flow cytometry, immunohistochemistry, western blotting, immunoprecipitation, ELISA, ELISpot, and multiplex assays.
<b>Expiration Date</b>	12 months from date of receipt.



Standard curve generated with Rat Anti-M...