


Product Datasheet

CSF3 Antibody (Biotin) (orb1250887)

Description	CSF3 Antibody (Biotin)
Species/Host	Rat
Reactivity	Human
Conjugation	Biotin
Tested Applications	Blocking, FC, ICC, IHC-Fr, IP
Immunogen	E. coli-expressed human G-CSF
Target	CSF3
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	G-CSF Antibody [BVD11-37G10] for use in ELISA and immunohistochemistry / immunocytochemistry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P09919
NCBI	P09919
Dilution Range	G-CSF Antibody [BVD11-37G10] for use in ELISA and immunohistochemistry / immunocytochemistry assays.
Expiration Date	12 months from date of receipt.



Standard curve generated with Rat Anti-H...