


Product Datasheet

IL3 Antibody (Biotin) (orb1250912)

Description	IL3 Antibody (Biotin)
Species/Host	Rat
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC
Immunogen	Yeast-expressed human IL-3
Target	IL3
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	IL-3 Antibody [BVD3-1F9] for use in flow cytometry, immunohistochemistry, ELISA, western blotting, and immunoprecipitation assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P08700
NCBI	P08700
Dilution Range	IL-3 Antibody [BVD3-1F9] for use in flow cytometry, immunohistochemistry, ELISA, western blotting, and immunoprecipitation assays.
Expiration Date	12 months from date of receipt.



Standard curve generated with Rat Anti-H...