


Product Datasheet

IFNG Antibody (Biotin) (orb1250917)

Description	IFNG Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human, Monkey, Primate
Conjugation	Biotin
Tested Applications	ELISA, FC, ICC, IHC-Fr, IHC-P, IP, WB
Immunogen	E. coli-expressed human IFN-gamma
Target	IFNG
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	IFN-gamma Antibody [B27] for use in flow cytometry, immunohistochemistry / immunocytochemistry, ELISA, and immunoprecipitation assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P01579
NCBI	P01579
Dilution Range	IFN-gamma Antibody [B27] for use in flow cytometry, immunohistochemistry / immunocytochemistry, ELISA, and immunoprecipitation assays.
Expiration Date	12 months from date of receipt.



Standard curve generated with Mouse Anti...