


Product Datasheet

CSF2 Antibody (Biotin) (orb1250920)

Description	CSF2 Antibody (Biotin)
Species/Host	Rat
Reactivity	Human, Primate
Conjugation	Biotin
Tested Applications	ELISA, FC, ICC, IHC-Fr, IHC-P, IP, WB
Immunogen	E. coli-expressed human GM-CSF
Target	CSF2
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	GM-CSF Antibody [BVD2-21C11] for use in flow cytometry, immunohistochemistry / immunocytochemistry, ELISA, ELISpot, immunoprecipitation, and multiplex assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P04141
NCBI	P04141
Dilution Range	GM-CSF Antibody [BVD2-21C11] for use in flow cytometry, immunohistochemistry / immunocytochemistry, ELISA, ELISpot, immunoprecipitation, and multiplex assays.
Expiration Date	12 months from date of receipt.



Standard curve generated with Rat Anti-H...