

## Product Datasheet

### IL10 Antibody (orb1250925)

<b>Description</b>	IL10 Antibody
<b>Species/Host</b>	Rat
<b>Reactivity</b>	Human, Primate, Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC
<b>Immunogen</b>	COS-expressed human IL-10
<b>Target</b>	IL10
<b>Preservatives</b>	Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	0.5 mg/mL
<b>Storage</b>	Store vial at 2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	IL-10 Antibody [JES3-12G8] for use in immunohistochemistry, ELISA, ELISpot, and western blotting assays.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P22301</b>
<b>NCBI</b>	<b>P22301</b>
<b>Dilution Range</b>	IL-10 Antibody [JES3-12G8] for use in immunohistochemistry, ELISA, ELISpot, and western blotting assays.
<b>Expiration Date</b>	12 months from date of receipt.