


Product Datasheet

B2M Antibody (FITC) (orb1250976)

Description	B2M Antibody (FITC)
Species/Host	Mouse
Reactivity	Aves, Gallus
Conjugation	FITC
Tested Applications	ELISA, FC, IP
Immunogen	Affinity purified chicken MHC-I molecules
Target	B2M
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	beta2-Microglobulin Antibody [F21-21] for use in flow cytometry, immunocytochemistry, and western blotting assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P21611
NCBI	P21611
Dilution Range	beta2-Microglobulin Antibody [F21-21] for use in flow cytometry, immunocytochemistry, and western blotting assays.
Expiration Date	12 months from date of receipt.



Chicken peripheral blood mononuclear cel...