


Product Datasheet

MHC Class I Antibody [F21-2] (Biotin) (orb1250981)

Description	MHC Class I Antibody [F21-2] (Biotin)
Species/Host	Mouse
Reactivity	Aves, Gallus
Conjugation	Biotin
Tested Applications	FC
Preservatives	Supplied in PBS/NaN ₃
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	MHC Class I Antibody [F21-2] for use in flow cytometry, immunohistochemistry, electron microscopy, western blotting, and immunoprecipitation assays.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	MHC Class I Antibody [F21-2] for use in flow cytometry, immunohistochemistry, electron microscopy, western blotting, and immunoprecipitation assays.
Expiration Date	12 months from date of receipt.



Chicken peripheral blood mononuclear cel...