


Product Datasheet

IgY Antibody [G-1] (FITC) (orb1250991)

| | |
|----------------------------|-----------------------------------------------------------------------------------------------------|
| Description | IgY Antibody [G-1] (FITC) |
| Species/Host | Mouse |
| Reactivity | Aves, Gallus |
| Conjugation | FITC |
| Tested Applications | FC, IHC-Fr |
| Preservatives | Supplied in PBS/NaN3 |
| Form/Appearance | Liquid |
| Concentration | 0.5 mg/mL |
| Storage | Store vial at 2-8°C |
| Note | For research use only |
| Application notes | IgY Antibody [G-1] for use in flow cytometry and immunohistochemistry / immunocytochemistry assays. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Dilution Range | IgY Antibody [G-1] for use in flow cytometry and immunohistochemistry / immunocytochemistry assays. |
| Expiration Date | 12 months from date of receipt. |



FLISA plate was coated with purified chi...