

Product Datasheet

S-Tag Antibody [SBSTAGa] (FITC) (orb1251045)

Description	S-Tag Antibody [SBSTAGa] (FITC)
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC, IHC-Fr, IP
Immunogen	Recombinant S-tag protein
Preservatives	Supplied in PBS/NaN ₃
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	S-Tag Antibody [SBSTAGa] for use in western blotting and other fluorescent assays.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	S-Tag Antibody [SBSTAGa] for use in western blotting and other fluorescent assays.
Expiration Date	12 months from date of receipt.