

Product Datasheet

S-Tag Antibody [SBSTAGa] (orb1251046)

Description	S-Tag Antibody [SBSTAGa]
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Recombinant S-tag protein
Preservatives	Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	S-Tag Antibody [SBSTAGa] for use in ELISA and western blotting assays.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	S-Tag Antibody [SBSTAGa] for use in ELISA and western blotting assays.
Expiration Date	12 months from date of receipt.