


Product Datasheet

Icam2 Antibody (FITC) (orb1251079)

Description	Icam2 Antibody (FITC)
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	COS cells transfected with mouse ICAM-2 cDNA and BW5147 AKR/J mouse thymoma cell line
Target	Icam2
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD102 Antibody [3C4] for use in flow cytometry and immunohistochemistry / immunocytochemistry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P35330
NCBI	P35330
Dilution Range	CD102 Antibody [3C4] for use in flow cytometry and immunohistochemistry / immunocytochemistry assays.
Expiration Date	12 months from date of receipt.



BALB/c mouse splenocytes were stained wi...