

Product Datasheet

LPPR3 antibody (orb125108)

Description	Goat Polyclonal to LPPR3. NP_079164.1, PMID: 12730698 and GeneID: 79948.
Species/Host	Goat
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Target	LPPR3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments in Human Brain (Cerebellum, Cerebral Cortex and Hippocampus) lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results
Protein Sequence	TQRGHDSVYQQN
Clonality	Polyclonal
MW	78.9
Entrez	79948
Dilution Range	ELISA: 1:32000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.